




THD REUSABLE NEEDLE HOLDER



Use

The THD Reusable Needle Holder is used for the treatment of the Haemorrhoidal disease, using the Transanal Haemorrhoidal Dearterialisation technique.

General Description

THD Reusable Needle Holder , is a reusable, steel needle holder with tungsten carbide inserts, and includes laser etched line at the tip of the holder to indicate proper needle placement.

Sterilization

Autoclave: Standard Surgical tool cycle (134°C)

Dimensions	
Length	200 mm

THD SUTURE		
Qty	P/N	Description
1	880003	THD Reusable Needle Holder